Send Orders for Reprints to reprints@benthamscience.ae

The Open Biotechnology Journal, 2017, 11, 89-93



Content list available at: www.benthamopen.com/TOBIOTJ/

DOI: 10.2174/1874070701711010089

RESEARCH ARTICLE

An Easy Method for Screening and Detection of Laccase Activity

Albino A. Dias^{*}, António J.S. Matos, Irene Fraga, Ana Sampaio and Rui M.F. Bezerra

Centre for the Research and Technology of Agro-Environment and Biological Sciences - CITAB, Universidade de Trásos-Montes e Alto Douro – UTAD, 5000-801 Vila Real, Portugal

Received: March 14, 2017

Revised: July 21, 2017

Accepted: August 13, 2017

Abstract: **Objective:**

An instrument-free assay was developed for simultaneous detection of laccase activity in a large number of samples as diverse as screening of laccase-producing microbial cultures or chromatographic fractions.

Method:

Dried paper discs previously impregnated with 2,2'-azino-bis(3-ethylbenzothiazoline-6-sulfonic acid) were placed on a flat-bottom microplate (a simple way to avoid misidentification) and loaded with an aliquot from each sample.

Results:

Discs corresponding to samples containing laccase activity become green-bluish colored within first ten minutes of reaction, allowing direct detection through simple naked-eye inspection.

Conclusion:

As an example, this easy process was applied to the laccase purification in order to eliminate chromatographic fractions that did not contain laccase activity, thus reducing the number of spectrophotometric assays.

Keywords: ABTS, Laccase, White-rot fungi, Zymogram, Laccase purification, Laccase activity.

1. INTRODUCTION

Laccases (benzenediol: oxygen oxidoreductase, EC 1.10.3.2) and laccase-like metalloxidases [1] are glycoenzymes that belong to the family of multicopper oxidases. They catalyze non-specific oxidation of natural or synthetic substrates such as phenolic or non-phenolic aromatic compounds, amines and inorganic ions or molecules, with a concomitant four-electron reduction of molecular oxygen to water. In recent years, these enzymes have received increasing attention due to their potential for wide-ranging biotechnological applications [2, 3] including environmental pollution abatement [4, 5]. As a consequence, a growing number of laccases from diverse organisms (filamentous fungi, yeasts, plants, bacteria, and even insects) but mainly from species belonging to the ecological group of white-rot basidiomycetes have been characterized at biochemical and molecular levels [2, 6].

As a rule of thumb, a process of isolation and purification or at least partially enzyme purification [7] should be done before proceeding with detailed structure-function characterization including determination of kinetic parameters. More than one chromatographic step is usually performed during purification procedures in order to obtain an enzymeenriched preparation free of other proteins or contaminant compounds such as inhibitors. In general, chromatographic processes yield a large number of fractions that have to be analyzed for the presence of enzyme activity, usually by

89



Address correspondence to this author at the Centre for the Research and Technology of Agro-Environment and Biological Sciences - CITAB, Universidade de Trás-os-Montes e Alto Douro - UTAD, 5000-801 Vila Real, Portugal, Tel: +351 259350250; E-mail: jdias@utad.pt

spectrophotometry. The method described here is a simple screening procedure that enables rapid detection of laccasecontaining samples. In this work, its usefulness is illustrated, but not limited, to detection of laccase activity in chromatographic fractions during purification of laccase produced by the white-rot fungus *Ganoderma resinaceum*.

2. MATERIALS AND METHODS

2.1. Production and Purification of Fungal Laccase

The fungal strain *G. resinaceum* Ga-20, from the culture collection at Laboratory of Agro-Environmental Technologies of UTAD, was grown in Potato Dextrose Agar (PDA) plates and incubated at 28°C for six days. As described earlier [8], Erlenmeyer flasks containing 100 ml of liquid medium supplemented with 5% (v/v) olive mill wastewater (whose main characteristics were as follows: 0.8 g phenols l^{-1} , 12.6 g total solids l^{-1} , chemical oxygen demand 53.5 g O₂ l^{-1}) were inoculated with two 1 cm² agar plugs removed from fungal PDA plates and incubated at 28°C and 120 rpm.

Culture medium was harvested at 12^{th} day of incubation, freeze-thawed, filtered (Whatman GF/A) and simultaneously equilibrated with 15 mM Tris-HCl pH 7.0 (buffer A) and concentrated by ultrafiltration performed under nitrogen at a positive pressure of 3.5 bar in an Amicon stirred cells unit fitted with a 10 kDa cut-off membrane. Subsequent anion exchange chromatography of concentrated crude laccase was accomplished in a DEAE-Sepharose FF column (10 x 150 mm) connected to a Gilson system (two peristaltic pumps and dynamic mixer, UV/Vis detector, and fraction collector). After baseline stabilization, bound proteins were eluted at a flow rate of 1.6 ml min⁻¹ with the following gradient of buffer B (15 mM Tris-HCl pH 7.0; 1 M NaCl): t = 0 min, %B = 0; t = 20 min, %B = 0; t = 54 min, %B = 50; t = 58 min, %B = 100; t = 60 min, %B = 100; t = 61 min, %B = 0. The fraction collector was set to change fraction tubes at 2-min intervals, which means a total of 30 fractions.

2.2. Easy Assay for Rapid Screening of Laccase Activity

2.2.1. Materials and Its Preparation

Standard 6 mm blank paper discs (BBL; Becton Dickinson & Company) were soaked for 2 min in a beaker with a solution of 2 mM 2,2'-azino-bis(3-ethylbenzothiazoline-6-sulfonic acid) (ABTS) prepared in 100 mM phosphate-citrate buffer pH 4.0. Subsequently, ABTS-wet discs were removed with the aid of a simple colander and dried in an oven at 60°C. After cooling, dried ABTS-impregnated discs should be stored in a sealable glass or plastic flask in the refrigerator until use.

2.2.2. Methodology

In this work, 10 dried ABTS-impregnated discs per each row were placed into an empty standard flat-bottom 96well microplate. The screening of laccase activity was started by the addition of 10 µl aliquots from each sample (in this work, chromatography fractions) to discs. Within the first 10 min of incubation at 30°C, discs loaded with aliquots from chromatography fractions containing laccase activity will present a green-bluish color development, easily visible by simple naked eye inspection.

2.3. Determination of Laccase Activity

Spectrophotometric determination of laccase activity was carried out at 25°C by measuring the oxidation of 2 mM ABTS in 100 mM phosphate-citrate buffer pH 4.0. The formation of cation radical ABTS⁺⁺ was followed by absorbance measured at 420 nm ($\varepsilon_{420} = 36 \text{ mM}^{-1} \text{ cm}^{-1}$) using an UV-Vis spectrophotometer [9]. One enzyme unit (U) is defined as the amount of enzyme which oxidizes 1 µmol of substrate per minute in the assay conditions.

3. RESULTS AND DISCUSSION

The laccase activity in culture media fermented by the white-rot fungus used in this work (*G. resinaceum*) or other laccase-producing white-rot fungal strains, such as *Trametes versicolor*, *Phlebia rufa* and *Fomes fomentarius*, are routinely detected in our laboratory using the proposed easy assay for rapid screening of laccase activity. Thus, the use of *G. resinaceum* in this study "*per se*" is not important or relevant to this work, since it was only used as a laccase-producing strain.

In this work, taking into account the chromatographic runtime and the fact that the fraction collector was set to

Fast Screening and Detection of Laccase Activity

collect and change the tubes at every two minutes intervals we have about 30 fractions per each chromatographic run that need to be analyzed for the presence of laccase activity. In order to cut down the number of spectrophotometric measurements of laccase activity, a simple zymogram assay was developed which relies on an equipment-free combination of two products readily available in most laboratories: blank paper discs (previously impregnated with ABTS) and microplates.

The upper panel of Fig. (1) shows the results of the proposed zymogram assay (color development of ABTSimpregnated paper discs) and the lower panel of Fig. (1) shows the plot of laccase activity of selected corresponding chromatographic fractions measured by the spectrophotometric assay (fractions number 2, 7, 10, 27 and 29 were not measured by spectrophotometry). As can be seen, most of the laccase activity occurs between fractions 19 - 21 as detected by the proposed zymogram assay. In fact, the peak of laccase activity was eluted with a saline gradient at about 0.28 M NaCl and was recovered in fraction number 20 (elution time between 38 - 40 min). Since this result matches the spectrophotometric determinations, proves that the zymogram assay is a useful, fast and efficient screening procedure.




Fig. (1). Screening of laccase activity in chromatography fractions recovered from an anion exchange column (DEAE-Sepharose FF) using a microplate containing dried discs impregnated with ABTS. Upper panel: row A: fractions 1 - 10; row B: fractions 11 - 20; row C: fractions 21 - 30. Lower panel: corresponding chromatography graph showing the elution profile of total proteins (—) and spectrophotometric determination of laccase activity (•) in selected fractions (row A fractions 2, 7and 10 as well as row C fractions 27 and 29 were not determined by spectrophotometry).

Diverse techniques have been used for screening of extracellular lignocellulolytic enzymes and for detection of laccase in purified samples [10 - 12]. The proposed method described here, present the following main advantages: i) blank paper discs commonly used for compound sensitivity tests were chosen due to its appropriated dimensions such as thickness and friendly handling; ii) paper discs were impregnated with a diluted ABTS solution (2 mM) but it should

be buffered at an acidic pH (*e.g.* 4.0); iii) the use of a microplate as a support for the ABTS-impregnated discs, also is a simple way to match each disc to a given fraction, thus avoiding sample misidentification.

The use of a buffered substrate solution is a critical point that needs particular emphasis because using paper discs impregnated with non-buffered ABTS solution can lead to false negative results. This is due to a monotonic decrease of laccase activity as the reaction pH value increases, typically observed in the pH-activity profile using ABTS as substrate [13]. In fact, a pH between 3.0- 4.0 is the optimal reaction value for the oxidation of ABTS catalyzed by most fungal and bacterial laccases [13 - 15]. Finally, dried paper discs impregnated with substrate can be stored in a refrigerator for long time, ready to be used when needed (according to our own experience no activity loss was observed after six months storage).

CONCLUSION

A new zymogram assay was applied to the detection of laccase activity. As a concluding remark, the method proposed here is a simple and economic process that will greatly cut down on tedious and unnecessary spectrophotometric assays of laccase activity during the screening of large number of samples as is the case of purification steps or detection of laccase activity in fermentation broths.

ETHICS APPROVAL AND CONSENT TO PARTICIPATE

Not applicable.

HUMAN AND ANIMAL RIGHTS

No Animals/Humans were used for studies that are base of this research.

CONSENT FOR PUBLICATION

Not applicable.

CONFLICT OF INTEREST

The authors declare no conflict of interest, financial or otherwise.

ACKNOWLEDGEMENTS

Declared none.

REFERENCES

- Martins LO, Durão P, Brissos V, Lindley PF. Laccases of prokaryotic origin: enzymes at the interface of protein science and protein technology. Cell Mol Life Sci 2015; 72(5): 911-22.
 [http://dx.doi.org/10.1007/s00018-014-1822-x] [PMID: 25572294]
- [2] Giardina P, Faraco V, Pezzella C, Piscitelli A, Vanhulle S, Sannia G. Laccases: a never-ending story. Cell Mol Life Sci 2010; 67(3): 369-85. [http://dx.doi.org/10.1007/s00018-009-0169-1] [PMID: 19844659]
- Zhao L, Yin W, Wang L. High-yield laccase-producing strains constructed by protoplast fusion between bacterium and fungus. Open Biotechnol J 2015; 9: 221-4. [http://dx.doi.org/10.2174/1874070701509010221]
- Majeau JA, Brar SK, Tyagi RD. Laccases for removal of recalcitrant and emerging pollutants. Bioresour Technol 2010; 101(7): 2331-50. [http://dx.doi.org/10.1016/j.biortech.2009.10.087] [PMID: 19948398]
- [5] Strong PJ, Claus H. Laccase: a review of its past and its future in bioremediation. Crit Rev Environ Sci Technol 2011; 41: 373-434. [http://dx.doi.org/10.1080/10643380902945706]
- [6] Singh A. D, Kumar S. R. Ligninolytic fungal laccases and their biotechnological applications. Appl Biochem Biotechnol 2010; 160(6): 1760-88.

[http://dx.doi.org/10.1007/s12010-009-8676-y] [PMID: 19513857]

- [7] Singha S, Panda T. Optimization of laccase fermentation and evaluation of kinetic and thermodynamic parameters of a partially purified laccase produced by *Daedalea flavida*. Prep Biochem Biotechnol 2015; 45(4): 307-35. [http://dx.doi.org/10.1080/10826068.2014.887581] [PMID: 24547974]
- [8] Matos AJ, Bezerra RM, Dias AA. Screening of fungal isolates and properties of *Ganoderma applanatum* intended for olive mill wastewater decolourization and dephenolization. Lett Appl Microbiol 2007; 45(3): 270-5. [http://dx.doi.org/10.1111/j.1472-765X.2007.02181.x] [PMID: 17718838]

- [9] Dias AA, Bezerra RM, Pereira AN. Activity and elution profile of laccase during biological decolorization and dephenolization of olive mill wastewater. Bioresour Technol 2004; 92(1): 7-13.
 [http://dx.doi.org/10.1016/j.biortech.2003.08.006] [PMID: 14643980]
- [10] Sundman V, Näse L. A simple test for direct visualization of biological lignin degradation. Paper Timber 1971; 53: 67-71.
- [11] Marr CD. Laccase and tyrosinase oxidation of spot test reagents. Mycotaxon 1979; 9: 244-76.
- [12] Dekker RF, Ling KY, Barbosa AM. A simple method for monitoring chromatography column eluates for laccase activity during enzyme purification. Biotechnol Lett 2000; 22: 105-8. [http://dx.doi.org/10.1023/A:1005680027798]
- [13] Xu F. Effects of redox potential and hydroxide inhibition on the pH activity profile of fungal laccases. J Biol Chem 1997; 272(2): 924-8. [http://dx.doi.org/10.1074/jbc.272.2.924] [PMID: 8995383]
- [14] Leitner C, Hess J, Galhaup C, et al. Purification and characterization of a laccase from the white-rot fungus Trametes multicolor. Appl Biochem Biotechnol 2002; 98-100: 497-507.
 [http://dx.doi.org/10.1385/ABAB:98-100:1-9:497] [PMID: 12018276]
- [15] Kaushik G, Thakur IS. Purification, characterization and usage of thermotolerant laccase from *Bacillus sp.* for biodegradation of synthetic dyes. Appl Biochem Microbiol 2013; 49: 352-9. [http://dx.doi.org/10.1134/S0003683813040169]

© 2017 Dias et al.

This is an open access article distributed under the terms of the Creative Commons Attribution 4.0 International Public License (CC-BY 4.0), a copy of which is available at: (https://creativecommons.org/licenses/by/4.0/legalcode). This license permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited.